ORDRE DE PRELEVEMENT DES TUBES

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ORDER OF BLOOD COLLECTION TUBES

Please respect the order of the additives in the tubes, it could affect results.

Reminder ²

Blood tubes filling order

Venipuncture	Invert	Capillary puncture
 Blood culture (Aerobic bottle first, then anaerobic)	5 to 10 times	 Capillary gaz
 Sodium citrate	4 times	 EDTA
 With/without clot activator, with or with/without SST	5 to 10 times	 Sodium or lithium heparin and other anticoagulants
 Sodium or lithium heparin and other anticoagulants	8 to 10 times	 With/without clot activator, with or with/without SST
 EDTA	8 to 10 times	Slowly mix blood with anticoagulant
 EDTA (rare metals)	8 to 10 times	
 Potassium oxalate/Sodium fluoride (glycolysis inhibitor)	5 to 10 times	
 Sodium citrate (3,8% for sedimentation rate Westergren)	5 to 10 times	

FOR OPTIMAL RESULTS OF ADDITIVE VERSUS BLOOD
INVERT THE TUBE DELICATELY UP AND DOWN THE NUMBER OF TIMES INDICATED ABOVE

IMPORTANT

- + Identify the specimen in **FRONT** of the patient
- + It is mandatory to double check the identification, for example :
 - + RAMQ **AND** first and last name of the patient.
 - + DOB **AND** first and last name of the patient.
- + **Write on the tubes & requisition, the time, the date and initials of the nurse who performed the specimen collection.**
- + Safely, dispose the contaminated material used for specimen collection.
- + Place the specimen in the "Biohazard" bag and seal it properly.
- + Place the requisition in the back pocket of the bag.

POSSIBLE CAUSE OF HEMOLYSIS

- + Gauge needle too small: Ideally use a size between **19G to 23G**.
- + Tourniquet kept in place more than a minute: Hemoconcentration, K, ALT, CK, proteins, LDH, Bili increased by 8 to 10%.
- + Repeated opening and closing of the hand: Hemoconcentration K, ALT, CK, proteins, LDH, Bili increased by 8 to 10%.
- + Tube shaken vigorously after specimen collection.
- + Sample taken from a hematoma site.
- + Presence of excessive alcohol on the puncture site.

OTHERS

- + Avoid using fine butterfly needles: the dead volume of the tubular falsifies the blood/anticoagulant ratio.
- + Favor the use of a 21G needle to a fine butterfly needle.
- + **Draw a discharge tube for every analysis, when a butterfly needle is used (in order to create a vacuum in the tubular).**
- + Never transfer a sample from one tube to another after blood collection.
- + **Let the sample stand upright for 30 minutes to allow blood to clot (gel/SST)**

Centrifuge according to the applicable recommendations

FAILURE TO COMPLY WITH THIS PROCEDURE MAY AFFECT THE QUALITY OF RESULTS

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